



M I D S I Z E S C O O T E R

Seat Color
Black Vinyl



Body Colors
Blue



Ride

GOLDEN COMPANION

3 wheel (Midsize)

The all-new Companion scooter provides great style with high performance and user-friendly disassembly. New features include an adjustable LED headlight; new ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion offers great value in a full-time portable scooter.



Rise.

Roll.

Rest.

Ride.



GOLDEN
technologies

YOUR LIFE IN MOTION


Companion 3 wheel Specifications

Model	GC240
Medicare Code	K0806
Weight Capacity	300 lbs
Drive Wheels	Rear
Maximum Speed	4.5 mph
Operating Range ¹	15.5 miles
Ground Clearance	3.75"
Turning Radius	39"
Type Batteries	2-U1
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
Dimensions: ³	
Length	42.5"
Width	23.5"
Height (ground to top of back of seat)	36.5" - 39.5"
Ground to Top of Seat	23" - 26"
Ground to Top of Deck	6.25"
Deck to Top of Seat	17" - 20"
Number of Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	30"
Weight of Unit: (assembled)	169 lbs
Front Half	46 lbs
Rear Half	48 lbs
Seat w/Arms	28 lbs
Batteries	22 lbs ea
Tires:	Solid, low profile
Front	9"
Rear	9"
Rear Anti-Tippers	3"
Standard Seat:	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" X 16"
Color	Black Vinyl
Optional Seats: ³	Stadium Style
Back Height - no headrest	16"
Width x Depth	20"x 18"
Color	Black
Seat Options:	
Pneumatic Rise Seat Option	Yes
Electric Rise Seat Option	Yes
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years

¹. Battery range will vary due to rider weight, drive surface and terrain.

². Electronics warranty excludes batteries.

³. There is a +1"/-1" tolerance on all dimensions.



Angle adjustable LED headlight creates a wide path of bright light



EZ Reach tiller adjustment




Patented infinite adjustable tiller



Quiet & maintenance free motor

Standard Features on Every Model:

- Under carpet troubleshooting guide
- Delta tiller with auto lock mechanism
- EZ Reach tiller adjustment lever
- LiquiCell Comfort Equipped seat
- Seat slide adjustment lever
- Adjustable LED headlight



Blue Red



Disassembles easily!

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

