



LIQUICELL®
COMFORT EQUIPPED

Shown with optional
power elevating seat

F U L L S I Z E S C O O T E R

Seat Color
Black Vinyl



Body Colors
Blue



Ride

**GOLDEN
COMPANION**

3 wheel

The Companion three-wheel scooter offers more leg room and foot room with a wider front deck and a stylish new design. Features include easy disassembly; an adjustable LED headlight; ergonomic control panel design with backlit battery gauge display; and wraparound Delta tiller. The Companion is the ultimate value in a full-time portable scooter.



Rise.

Roll.

Rest.

Ride.



GOLDEN
technologies


YOUR LIFE IN MOTION

Companion 3 wheel Specifications


Model	GC340
Weight Capacity	350 lbs
Drive Wheels	Rear
Maximum Speed	4.5 mph
Operating Range ¹	16.5 miles
Ground Clearance to center deck	4.5"
Ground Clearance to transaxle/motor	2.5"
Turning Radius	40.5"
Type Batteries	2-U1
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	90 Amp Dynamic R Series
Dimensions:	
Length	45.5"
Width	23.5"
Height (ground to top of back of seat)	36.5" - 39.5"
Ground to Top of Seat	23" - 26"
Ground to Top of Deck	6.25"
Deck to Top of Seat	17" - 20"
Number of Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	34.5"
Weight of Unit: (assembled)	173 lbs
Front Half	48 lbs
Rear Half	48 lbs
Seat w/Arms	28 lbs
Batteries	22 lbs ea.
Tires:	Solid/Low Profile
Front	9"
Rear	10"
Rear Anti-Tippers	3"
Standard Seat:	Stadium Style
Slide Seat	Standard
Back Height - no headrest	16"
Width x Depth	18" X 16"
Color	Black Vinyl
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years

¹. Battery range will vary due to rider weight, drive surface and terrain.

². Electronics warranty excludes batteries.



Angle adjustable LED headlight creates a wide path of bright light



Quiet & maintenance free motor



EZ Reach tiller adjustment



Patented infinite adjustable tiller



Disassembles easily!

Standard Features on Every Model:

- Under carpet troubleshooting guide
- Delta tiller with auto lock mechanism
- EZ Reach tiller adjustment lever
- LiquiCell Comfort Equipped seat
- Seat slide adjustment lever
- Adjustable LED headlight

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.