



Golden Scooters

SUPERIOR DESIGNS | FIRST CLASS LUXURY



COMPACT
PORTABLE



LARGE
PORTABLE



FULL-SIZE
LUXURY



BUZZaround®

*Buzzaround
Lite
GB106*



*Buzzaround XL
GB116*

*Buzzaround XL
GB146*

COMPACT PORTABLE

The Buzzaround Series is perfect for those individuals seeking a fast and convenient way to travel. Do you have a long day of walking planned? Not a problem for the Buzzarounds, simply assemble/disassemble and get ready for a day full of cruising light.

BUZZAROUND XL MODELS

- 300 lb. weight capacity
- Delta Tiller standard
- Disassembles into 5 lightweight pieces
- Battery charging port conveniently located below the control panel
- Top speed 4.5 mph
- Wireless and easy disassembly!
- Comfortable 17"x16" seat standard
- Two sets of interchangeable shroud panels **INCLUDED FREE:** Red and Blue
- Under carpet troubleshooting guide included

BUZZAROUND LITE

- 250 lb. weight capacity
- Delta Tiller standard
- Disassembles into 5 lightweight pieces
- Top speed 4 mph
- Three sets of interchangeable shroud panels **INCLUDED FREE:** Red, Blue and Sahara



*Companion
GC340*

GOLDEN COMPANION



FULL-SIZE LUXURY

The Golden Companion scooters are your ticket to ride first class every day! Allow the luxurious, full-size Companions to enhance your life and take you to the places you've always had trouble getting to. Shop in the mall, attend sporting events and enjoy the great outdoors with the unmatched comfort of your favorite new Companion! All Companion models feature:

- All new design for easy disassembly
- Wraparound Delta Tiller standard
- LiquiCell™ comfort equipped seat
- Adjustable LED headlight
- EZ Reach tiller adjustment lever
- Auto lock tiller
- Easy to read control panel
- Easy access charging port
- Backlit battery gauge
- Black, non-marking, low profile tires on sturdy steel rims
- Comfortable, fully adjustable padded armrests
- Under-carpet troubleshooting guide
- Removable shopping basket with handle
- Four mobility transportation tie downs



Companion GC440

LiteRider[®] GOLDEN



LiteRider GL110



LiteRider GL140

LARGE PORTABLE

The unique design of the LiteRider scooters gets you wherever you need to go in a style all their own. The only scooters to come with two cargo baskets standard, the LiteRiders are easy to disassemble when needed, and maneuver smoothly indoors and out. Both LiteRider models feature:

- 300 lb. weight capacity
- Wraparound Delta Tiller standard
- Adjustable tiller standard
- Easy to read control panel with backlit battery gauge
- Two storage baskets standard
- Comfortable 17"x16" stadium style seat
- Adjustable armrests
- 18Ah batteries
- Completely wireless disassembly
- 9" front and rear tires
- Low-profile, black tires on high-strength steel rims with stylish wheel covers



Companion GC240



golden **Avenger**

FULL-SIZE RUGGED LUXURY

When you're in need of a rugged, heavy duty scooter to tackle the challenging terrain of the great outdoors, look no further than the Avenger. Featuring oversized tires with sporty mag wheels, a full-size captain's seat, and a top speed of 7mph, the Avenger will take you through life's biggest adventures. The Avenger model features:

- Adjustable tiller for riding comfort
- Deluxe, reclining high back captain's seat
- Large, 13" front and rear tires
- Adjustable front headlight and four-way flashers
- Large front cargo basket standard
- 18-mile operating range

AVAILABLE ACCESSORIES



Comfort Ride*



Safety Flag



Crutch Holder



Rear Basket



Oxygen Holder



Walker Holder



Cup Holder



Pack'n'Go



Power Elevating Seat*



Quad Cane Holder

*Only available on Companion models.

SPECIFICATIONS*

| Scoter Model | Length | Width | Tires | Front Tire | Rear Tire | Weight Capacity (lbs.) | Top Speed (mph) | Ground Clearance | Operating Range** (miles) | Turning Radius | Batteries (2) | Weight w/ Batteries | Front Half Weight | Rear Half Weight | Seat Weight w/ Arms | Heaviest Piece | Freewheel Mode | Electronic Speed Control | Charger |
|--------------|--------|-------|-------------------|------------|-----------|------------------------|-----------------|------------------|---------------------------|----------------|---------------------|---------------------|-------------------|------------------|---------------------|----------------|----------------|--------------------------|-----------|
| GB-106 | 36.5" | 21.5" | Solid Low Profile | 8" | 8" | 250 | 4 | 1.75" | 8 | 31" | 2-12Ah | 96 | 29 | 26 | 22 | 29 | Yes | Yes | Off-board |
| GB-116 | 40.5" | 21.5" | Solid Low Profile | 8" | 8" | 300 | 4.5 | 2" | 7/10 ¹ | 37" | 2-12 Ah/ 2-18 Ah | 102 | 33 | 25 | 22 | 33 | Yes | Yes | Off-board |
| GB-146 | 42" | 21.5" | Solid Low Profile | 8" | 8" | 300 | 4.5 | 2" | 7.5/11 ¹ | 48.5" | 2-12 Ah/ 2-18 Ah | 109 | 41 | 25 | 22 | 41 | Yes | Yes | Off-board |
| GL-110 | 41.5" | 22" | Solid Low Profile | 9" | 9" | 300 | 5 | 3.75" | 9.5 | 39" | 2-18Ah | 119 | 44 | 28 | 22 | 44 | Yes | Yes | Off-board |
| GL-140 | 43.25" | 22" | Solid Low Profile | 9" | 9" | 300 | 5 | 3.75" | 9.5 | 48" | 2-18Ah | 127 | 52 | 28 | 22 | 52 | Yes | Yes | Off-board |
| GC-240 | 42.5" | 23.5" | Solid Low Profile | 9" | 9" | 300 | 4.5 | 3.75" | 15.5 | 39" | 2-U1 | 166 | 46 | 48 | 28 | 48 | Yes | Yes | Off-board |
| GC-340 | 45.5" | 23.5" | Solid Low Profile | 9" | 10" | 350 | 4.5 | 4.5" | 16.5 | 40.5" | 2-U1 | 168 | 48 | 48 | 28 | 48 | Yes | Yes | Off-board |
| GC-440 | 47.5" | 23.5" | Solid Low Profile | 10" | 10" | 350 | 4.5 | 4.5" | 15 | 53.5" | 2-U1 | 183 | 63 | 48 | 28 | 63 | Yes | Yes | Off-board |
| GA-541 | 57.5" | 24.5" | Pneumatic | 13" | 13" | 500 | 7 | 4" | 18 | 67" | 2-NF22 | 265 | 68 | 68 | 55 | 68 | Yes | Yes | Off-board |

* ALL WEIGHT SHOWN IN LBS.

** Battery range will vary due to rider weight, drive surface and drive slope. Range was tested at maximum weight.

¹. Range tested at 200 lbs.

Brakes: All Golden Technologies scooters come standard with Dynamic/Electromechanical braking systems. **Rear wheel drive:** All Golden Technologies scooters come standard with a sealed transaxle, 24 volt motor.