

Supernal™ Recliner Beds



Head Tilt Position



Head Fully Elevated



Flat Fully Elevated



Wireless Hand Control



Fully Elevated

Features:

- Illuminated Wireless Hand Control
- Massage system
- Wall Hugging action
- DC powered motor and control box
- Emergency Battery Power Down
- Several mattress options
- Low Profile Floor guides
- Vertical Lift is straight up for 10 inches
- Optional half-length side rails available

Specifications:

Angles	Head: 65 Degrees Foot: 35 Degrees
Vertical Heights <i>(Deck to floor)</i>	Lowest: 10 inches Highest: 21 inches
Outer Dimension	Twin: XL 38"x79" Full: 52"x79" Queen: 59"x79"

Weight Capacity	Twin: 400lbs Full: 500lbs Queen: 600lbs
------------------------	---

Frame Weights	Recliner Plus	Hi-Low Module
Twin:	165lbs	112lbs
Full:	227lbs	136lbs
Queen:	232lbs	135lbs



Distributed By:
Golden Technologies of Canada
740 Huronia Road, Unit 5
Barrie, ON, L4N 6C6
Tel: (800) 265-9615
www.goldentechnologies.ca



Supernal™ Mattress Options

The Supernal Sleep System is a one-stop shop for all of your mattress needs. The two mattresses in this line, the Ascent 1 and the Soft Touch are the perfect matches for our Supernal Hi-Low beds, as well as an outstanding choice for replacing an existing adjustable bed. The Ascent 1 is constructed of state-of-the-art progressive laminated foams to provide a healthy skin environment. Designed with sleep technologies that promote excellent comfort and healthy bedding solutions. Manufactured in an FDA approved facility that has decades of manufacturing experience to ensure the highest quality. The Soft Touch mattress is made of cool, form-fitting memory foam that conforms to the shape of the person who is using it. Memory foam users credit their wholesome night's rest to the absolute comfort of the soft memory foam. Both mattresses come with border edge protection for added comfort and peace of mind.

Two mattress cover options are available the Natural Bamboo Quilted Cover and the Stretch Fabric Cover. The Natural Bamboo Quilted Cover provides a comfortable and fresh sleeping environment. Natural properties: anti-bacterial, deodorizing, and breathable. Comprised of Joma wool filling, which absorbs and causes moisture to evaporate, for maximum comfort. The Stretch Fabric Cover is a unique, waterproof cover that feels soft, just like fabric. It can be cleaned with household cleansers. It will protect your mattress for many years to come.

Specifications:

Outer Dimension Twin: XL 38"x79"
Full: 52"x79"
Queen: 59"x79"

Mattress Height 6 inches

Mattress Weights Twin: 43lbs
Full: 59lbs
Queen: 67lbs




Removable cover



Soft Touch memory foam



Natural Bamboo Quilted Cover



Stretch Fabric Cover



Distributed By:
Golden Technologies of Canada
740 Huronia Road, Unit 5
Barrie, ON, L4N 6C6
Tel: (800) 265-9615
www.goldentechnologies.ca

Supernal™
SLEEP SYSTEM