

V A L U E S E R I E S

Rise



Capri
PR-200



Autumn



Cornflower



Spring



Rosewood



Chestnut

Capri

The Capri's classic Colonial button-back design is enhanced by a slightly firmer lumbar support area. Its arms offer a wider, more rounded design for extra comfort and stability. A truly sensible lift/recliner, the Capri combines dependability with long-lasting comfort and style.



Rise.

Roll.

Ride.

GOLDEN
technologies

YOUR LIFE IN MOTION

PROUDLY

MADE IN AMERICA
with Domestic & Imported Materials

Capri



Capri

Model Number	PR-200
Size	Medium
Weight Capacity	325 Lb.
Recline Positions	2
Standard Fabrics	5
Fabric Options	NONE
Chaise Seat	NO
Coil Springs	NO
Back Type	Button
Head Rest Cover	YES
Arm Covers	YES
Battery Back Up	YES
Overall Width	31"
Overall Height	41"
Width Between Arms	20"
Floor to Top of Seat	18"
Seat Depth	19"
Seat to Top of Back	27"
Seat to Top of Arm	7"
Seat Angle Raised	16°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	35"
Back Angle Closed	110°
Back Angle Reclined	135°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	160°
Distance Required from Wall Reclined	11"
Extension Reclined from Back of Seat to End of Foot Rest	40"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	1 Yr.
Pro Rated Years 4 Thru 7	NO
Warranty on Mechanical Labor	6 Mos.
Warranty on Electrical Labor	6 Mos.
Removable Back Option Available	NO
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	YES
Extra Pocket Option Available	YES
Heat and Massage Option Available	NO
Foot Rest Extension Option Available	NO
Straight Lift Option Available	YES
User Height (suggested)	5'4" -5'10"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Battery Backup System

Battery Backup System

Batteries not included. Battery Backup System is designed to bring the chair from the reclined to the seated position in the event of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68
Listing # E112417
FNAE 80-214A

