



**small<sup>to</sup>TALL**  
WE BUILD THEM ALL

# LIFT CHAIR SIZING CHART

LIFT CHAIR SIZE	RECOMMENDED USER HEIGHT
PETITE/SMALL	UNDER 5'
SMALL	5' - 5'3"
SMALL/MEDIUM	5'2" - 5'6"
MEDIUM	5'4" - 5'10"
MEDIUM/LARGE	5'7" - 6'2"
LARGE	5'11" - 6'2"
TALL	6'3" - 6'8"

\*This chart is quick reference for sizing - for full sizing details see Lift Chair Sizing Guide