

New!



H E A V Y D U T Y P O V

- Seat Color
- Black **coming soon!**
- Body Colors
- Orange
 - Red

Golden
Patriot™
4 wheel

If you're looking for the ultimate outdoor driving experience, look no further than the Golden Patriot™. The Golden Patriot's rugged outdoor design features an off road package, four-wheel independent suspension, a full lighting package and a full sized captain's seat. With a 400 lb. weight capacity and 13" pneumatic tires, the Golden Patriot will give you the ride of your dreams.



Patriot 4 wheel Specifications


Model	GR575
Medicare Code	Not Coded
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	8 mph
Operating Range ¹	24 miles
Ground Clearance	6.75"
Turning Radius	76.5"
Type Batteries	2-75 AH
Colors	Red, Orange
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Mechanical Hand Brake	Yes
Charger	Off Board
Controller	140 Amp
Dimensions:	
Length	57"
Width	28"
Height (ground to top of back of seat)	46.5" - 48.5"
Height (ground to top of back of headrest)	54" - 59"
Ground to Top of Seat	27" - 29"
Ground to Top of Deck	8.75"
Deck to Top of Seat	18" - 20"
Number of Seat Height Positions	3
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight of Unit: (assembled)	308 lbs
Base	152 lbs
Seat w/Arms	50 lbs
Batteries	52 lbs ea
Tires:	Pneumatic
Front	13"
Rear	13"
Standard Seat:	Captain's Pillow Top High-Back
Slide Seat	Yes
Back Height - no headrest	20"
Back Height with headrest	27.5" - 30.5"
Width x Depth	20" x 18"
Color	Black (coming soon)
Warranty:	
Frame	5 Years
Drive Train	1 Year
Electronics ²	1 Year
Batteries	1 Year

¹ Battery range tested at 200 lbs., and will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.



New button covers, hand (parking) brake lock, LED battery meter.



Updated control panel on wraparound delta tiller.



New 13" pneumatic tires.



New ultrabright LED rear lights.