




Shown in Antler


Standard Fabrics:



Antler



Lakefront



Chestnut

This contemporary plush fabric comes in beautiful, rich colors. It features our exclusive SofTouch™ with Stain Defense coating.

Crossover fabrics are unavailable.

Elara
PL118



Comfort Zone™
by GOLDEN

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Elara is offered in two sizes with a back style that provides maximum comfort in a tasteful seam design. It features a manually adjustable three-position headrest, supportive chaise pad, and a USB port on the hand control. Available in luxurious new fabrics to match any décor. Enjoy the powerful energy and take time to better connect with yourself in your DeLuna recliner by Golden!

Standard Features:

- Manually Adjustable Three-Position Headrest
- Easy Grip Armrest
- Soft-Feel Hand Control
- USB Charging Port



DeLUNA Series PL118

Model Number	PL118-MSM	PL118-LAR
Size	Medium-Small	Large
Weight Capacity	375 LBS	375 LBS
Positioning	3	3
Comfort Zone(s)	1	1
Twilight Technology	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	YES
Back Type	Seam	Seam
Removable Back Type	Bracket	Bracket
Head Rest Cover	YES	YES
Arm Rest Covers	NO	NO
Overall Width	32"	34"
Overall Height	43"	46"
Width Between Arms	21"	23"
Floor to Top of Seat	18.5"	19.5"
Seat Depth	21"	21"
Seat to Top of Back	29"	31"
Seat to Top of Arm	7"	6"
Distance Required from Wall Reclined	19"	20"
Extension Reclined from Back of Seat to End of Foot Rest	42"	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Option	NO	NO
Heat & Massage Option	NO	NO
Articulating Headrest Option	Manual	Manual
Adjustable Lumbar Option	NO	NO
AutoFlex Option	NO	NO
Foot Rest Extension Option	NO	NO
Straight Lift Option	YES	YES
USB Charging Port	Hand Control	Hand Control
Battery Backup	YES	YES
Place of Origin	China	China



Elegant seamed back design and supple urethane fabric make this chair both pleasing to the eye and highly durable.

Shown in Chestnut



Shown in Antler

Headrest can be manually adjusted in three different positions to control comfort level and support.

SmartTek Diagnostic System STANDARD on all Golden models.

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.


Comfort Zone™ by GOLDEN

Comfort Zone 1 provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



LOWER BODY UPPER BODY HEADREST LUMBAR TWILIGHT TILT

WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies

SHT-PL118-011224