



Standard Fabrics:



Acorn



Graphite



Maple

This faux leather is a rich and modern alternative to natural leather that is extremely comfortable while offering the rugged durability of leather. Made from 100% polyester, the fabric is easy to care for ensuring functionality as well as elegance. This long-lasting fabric will complement a variety of decors while standing the test of time under heavy use.

Crossover fabrics are unavailable.

Dione
PL446



Comfort Zone[™]
by GOLDEN

The Dione is offered in both medium and large sizes and features a stylish, overstuffed seamed back and plush chaise pad. It comes standard with adjustable lumbar and headrest, programmable AutoDrive[™] hand control, and USB charging port for your phone or tablet. Also standard is our patented four Comfort Zone positioning system, which provides a stress-free and supportive reclining experience.

Standard Features:

- Patented MaxiComfort[®] Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable AutoDrive[™] Hand Control
- USB Charging Port
- Extra Pocket



DeLUNA Series PL446

Model Number	PL446-MED	PL446-LAR
Size	Medium	Large
Weight Capacity	375 LBS	375 LBS
Positioning	MaxiComfort	MaxiComfort
Comfort Zone(s)	4	4
Twilight Technology	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	YES	YES
Back Type	Seam	Seam
Removable Back Type	Bracket	Bracket
Head Rest Cover	NO	NO
Arm Rest Covers	NO	NO
Overall Width	34"	35"
Overall Height	43.5"	46"
Width Between Arms	22"	22.5"
Floor to Top of Seat	20.5"	22"
Seat Depth	21"	21.5"
Seat to Top of Back	28"	30"
Seat to Top of Arm	7"	7.5"
Distance Required from Wall Reclined	24"	26"
Extension Reclined from Back of Seat to End of Foot Rest	42.5"	44"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Option	Standard	Standard
Heat & Massage Option	NO	NO
Articulating Headrest Option	Standard	Standard
Adjustable Lumbar Option	Standard	Standard
AutoFlex Option	NO	NO
Foot Rest Extension Option	NO	NO
Straight Lift Option	YES	YES
USB Charging Port	Outside Arm	Outside Arm
Battery Backup	NO	NO
Place of Origin	China	China



Shown in Graphite



Convenient built-in USB charging port for your phone or tablet located on the outside arm.

SmartTek Diagnostic System STANDARD on all Golden models.
Designed and engineered in Old Forge, PA. Assembled in China using Golden components.




Comfort Zone™ by GOLDEN

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



LOWER BODY UPPER BODY HEADREST LUMBAR TWILIGHT TILT

WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374
Fax: 800-628-5165
www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



www.facebook.com/goldentechnologies