



Buzzaround GB220



Pack your dreams and go on any adventure with the incredibly portable Buzzaround CarryOn® HD scooter! The **ONLY** lightweight travel scooter designed with the strength and power to accommodate up to a 400 lb. weight capacity. With a sporty blacked-out style, front and rear LED lighting, and optional color shrouds, this unique scooter not only looks great but offers excellent maneuverability and range. The CarryOn HD features a brushless motor for efficient operation, large wheels, impressive ground clearance, full sized seat, and a top speed of 4.29 mph making any destination possible. It offers ample foot and leg space, an adjustable tiller, and a large ergonomic seat for ultimate comfort while on the go. Enjoy long-range coverage up to 22 miles on the airline-friendly 12.5 AH and spare 6.5 AH lithium-ion batteries, which can be removed and taken as a carry-on while traveling by air, sea, or train. The heaviest piece is only 52.5 lbs. (without seat and battery) and the height when folded is less than 14" for storage where space is limited. A kickstand with a rear basket allows for vertical storage and serves as spare battery holder. Carry on with life using the CarryOn HD scooter!



FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.

Standard Colors:



Velvet Red
(Installed)



Coastal Blue
(Included in Box)

Optional Colors:



Crisp White



Canary Yellow



Stone Grey



Orchid Pink



Army Green
(Matte)



Desert Tan

Sporty Fold-Flat Scooter!
Airline Friendly!
A&C USB Ports!

Standard Features:

- 400 lb. Weight Capacity
- Brushless Motor
- Full Sized Seat
- Front & Rear LED Lights
- Width Adjustable Armrests
- Textured Soft Touch Armrests
- Removable Front Storage Basket
- 12.5 AH Lithium-ion Battery with Travel Case & 6.5 AH Spare
- Kickstand with Rear Basket & Spare Battery Holder
- Easily Interchangeable Color Shrouds
- Cane Holder
- Cup Holder
- A & C USB Ports

**FREE Spare 6.5 AH
Lithium-ion Battery for Up
to 22 Mile Combined Range
with 12.5 AH Battery!
2-Year Warranty!**



Model	GB220
Medicare Code	Not Coded
Weight Capacity	400 lbs
Drive Wheels	Rear
Maximum Speed	4.29 mph
Operating Range ¹ 12.5 AH Battery	14.2 miles
Operating Range ¹ 6.5 AH Battery	7.4 miles
Total Combined Range ¹	21.6 miles
Ground Clearance Under Deck	4.63"
Ground Clearance to Transaxle/Motor	2.3"
Turning Radius	44.5"
Battery Type (Lithium Ion)	12.5 AH & 6.5 AH
Colors	Red, Blue, White, Yellow, Grey, Pink, Green & Tan
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	32 Amp
Dimensions	
Length	40"
Width	23.5"
Height (ground to top of back of seat)	39"
Maximum Height (folded without seat)	14"
Ground to Top of Seat	23.75"
Ground to Top of Deck	6.6"
Deck to Top of Seat	17"
Front Axle to Rear Axle	28.7"
Weight of Unit: (assembled)	77.5 lbs
Main Body	52.5 lbs
Seat w/Arms	20 lbs
Batteries	5 lbs each
Tires:	Solid/Low Profile
Front	7.25"
Rear	8.25"
Rear Anti-Tip Wheels	2"
Standard Seat	Stadium Style
Back Height - no headrest	15"
Width x Depth	16" x 18"
Color	Black Vinyl
Warranty:	
Frame	Lifetime
Drive Train	2 Years
Electronics ²	2 Years
Batteries	2 Years

¹. Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain, and battery type.
². Electronics warranty excludes batteries.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Convenient kickstand for small space storage

LIGHTWEIGHT HEAVY-DUTY PORTABLE SCOOTER




Folds easily for a compact profile of only 14" high (without seat)



Optional spare battery storage in rear basket



Front & rear LED lighting



Front carry handle makes travel easy



Conveniently located USB ports



LED battery indicator on the console



Phone: 800-624-6374
 Fax: 800-628-5165
www.goldentech.com

Find us on Facebook
www.facebook.com/goldentechnologies

LIT-SHT-GB220, Rev.A 061125