




### Standard Colors:



Silver



Red

## Avenger 4 wheel GA541

golden  
**Avenger**

The Golden Avenger four-wheel, heavy duty scooter provides an impressive blend of rugged outdoor maneuverability while giving a smooth, comfortable ride over a variety of nature's most challenging terrains. Featuring a 500 lb. weight capacity.

### Features:

- 500 lb weight capacity
- Speeds up to 7 mph
- Bright safety lights
- 22 NF batteries
- Up to 18 Mile Range
- 4-Way Flashers
- Great Warranty
- Oversized tires
- Oversized seat
- Removable basket




### GA-541

## Avenger 4 wheel Specifications

<b>Model</b>	<b>GA541</b>
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	7 mph
Operating Range <sup>1</sup>	29.6 miles
Ground Clearance to Center Deck	4"
Turning Radius	67"
Type Batteries	2-NF22
Colors	Silver, Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PGS 140 Amp
<b>Dimensions:</b>	
Length	57.5"
Width	24.5"
Height (ground to top of back of seat)	42.5" - 45.5"
Height (ground to top of back of headrest)	50" - 56"
Ground to Top of Seat	25" - 28"
Ground to Top of Deck	6.5"
Deck to Top of Seat	18.75" - 21.75"
Number of Seat Height Positions	4
Size of Increments	1"
Front Axle to Rear Axle	38.5"
<b>Weight of Unit: (assembled)</b>	<b>265 lbs</b>
Front Half	68 lbs
Rear Half	68 lbs
Base	N/A
Seat w/Arms	55 lbs
Batteries	37 lbs ea
<b>Tires:</b>	<b>Pneumatic</b>
Front	13"
Rear	13"
Casters	N/A
Rear Anti-Tippers	3"
<b>Standard Seat:</b>	<b>Captain's Seat</b>
Slide Seat	N/A
Back Height - no headrest	21" - 24"
Back Height with headrest	27.5" - 31.5"
Width x Depth	20" X 17"
Color	Black Vinyl
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	2 Years
Electronics <sup>2</sup>	2 Years
Batteries	1 Year



Superb Style



<sup>1</sup> Battery range will vary due to rider weight, drive surface and terrain. Foam-filled tires are optional. The operational mileage will be less if scooter is equipped with foam filled tires. Batteries not included.

<sup>2</sup> Electronics warranty excludes batteries.