



Standard Colors:



Red Velvet
(Installed)



Coastal Blue
(Included)

ALL NEW Sporty Blacked-Out Styling!

Great for Indoor & Outdoor Use!

Buzzaround 3-Wheel GB107 4-Wheel GB108



The new and improved Buzzaround LT scooter is the smallest, lightest, and most portable Buzzaround scooter to get you where you want to go in comfort and style! Offered in both 3-wheel and 4-wheel models, upgrades include a sleek and sporty blacked-out design, a lightweight padded stadium-style seat, contoured hand grips, and a longer adjustable tiller for the ultimate driving experience. The new Buzz LT features a color display on the console, clear-coated plastic color through shrouds (red installed and blue included in the box), scuff guard on the foot deck, and additional side reflectors for added safety. With a weight capacity of 300 lbs., it also has increased foot room, a convenient charging port on the tiller, and a rugged hard plastic basket. Enjoy its indoor maneuverability and smooth outdoor ride with 2.5" of ground clearance! **The GB108 4-W scooter offers our new MicroTurn™ technology with an impressive 41.75" turning radius to navigate even the tightest spaces with ease!**

Standard Features:

- 300 lb. weight capacity
- Lightweight, padded stadium-style seat
- Ergonomic hand grips
- Non-slip deck with scuff guard
- Clear-coated, durable color through shrouds
- Battery charging port on the tiller
- Side reflectors
- Disassembles quickly and easily for transport



FDA Class II Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.

microTURN
NEW ON GB108!



LIGHTWEIGHT PORTABLE SCOOTER

Model	GB107	GB108
Class I Medical Device	Yes	Yes
Weight Capacity	300 lbs	300 lbs
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range ¹	7 miles	7.4 miles
Ground Clearance Under Deck	2.5"	2.5"
Turning Radius	39"	41.75"
Type Batteries	2-12AH	2-12AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	50 Amp	50 Amp
Dimensions:		
Length	41.75"	44"
Width	19-21"	19-21"
Height (ground to top of seat back)	35" - 36.5"	35" - 36.5"
Ground to Top of Seat	21" - 22.5"	21" - 22.5"
Ground to Top of Deck	4"	4"
Deck to Top of Seat	15" - 17"	15" - 17"
Number of Seat Height Positions	3	3
Size of Increments	.75"	.75"
Front Axle to Rear Axle	33.5"	33.5"
Weight of Unit: (assembled)		
Front Half	28 lbs	37.5 lbs
Rear Half	20.5 lbs	20.5 lbs
Seat w/Arms	19 lbs	19 lbs
Batteries (pack with 2 batteries)	17.5 lbs	17.5 lbs
Tires:		
	Solid, Low Profile	Solid, Low Profile
Front	8"	8"
Rear	8"	8"
Rear Anti-Tippers	2"	2"
Standard Seat:		
Back Height - no headrest	15"	15"
Width x Depth	17.5"x15"	17.5"x15"
Color	Black Vinyl	Black Vinyl
Warranty:		
Frame	Lifetime	Lifetime
Drive Train	2 Years	2 Years
Electronics ²	2 Years	2 Years
Battery	1 Year	1 Year

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.
² Electronics warranty excludes batteries.

microTURN

NEW MicroTurn feature on GB108 offers an impressive 41.75" turning radius!



Sleek dash with comfort grip handles and conveniently located charging port on the tiller.



Padded, contoured seat with width-adjustable, flip up arm rests for customized comfort.



Quickly disassembles into five, lightweight pieces for easy transport!



Phone: 800-265-9615
www.goldentechnologies.ca

Find us on Facebook
www.facebook.com/goldentechcanada

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.