



GB149 Shown



Standard Colors:



Velvet Red
(Installed)



Coastal Blue
(Included in Box)

NEW Sporty Optional Colors:



Crisp White



Canary Yellow



Stone Grey



Orchid Pink



Army Green
(Matte)



Desert Tan

Sporty New LED Running Light on the 4-W LX!

Buzzaround 3 wheel GB119 4 wheel GB149



The Buzzaround LX is America's first full-time luxury travel scooter! Full-size comfort comes with the convenience of true portability on the Buzzaround LX, along with standard features that include a stylish quilted high back seat, an infinite adjustable tiller at your fingertips, and a USB charging port right under the dash board. The 4-Wheel model features a sporty new LED running light on the front for maximum safety and visibility. With a 375 lb. weight capacity, comfort suspension springs, 40.5" turning radius, and up to 18 mile range on two U1 batteries, the Buzzaround LX offers ample foot room to take you where you want to go in comfort and style!

Standard Features:

- Stylish, foldable, high back seat with armrests
- 375 lb. weight capacity
- Infinitely adjustable tiller with EZ reach adjustment lever
- Up to 18 mile range with two U1 batteries
- Ample foot room
- USB charging port
- Comfort spring suspension for a smoother ride
- Cell phone holder





LIGHTWEIGHT PORTABLE SCOOTER

Full-size comfort with the convenience of true portability!

Item	GB119	GB149
Class I Medical Device	Yes	Yes
Weight Capacity	375 lbs.	375 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	5 mph	5 mph
Operating Range ¹	18 miles	18 miles
Ground Clearance Under Deck	4"	4"
Turning Radius	40.5"	54"
Type Batteries	2-U1	2-U1
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	90 Amp	90 Amp

Dimensions		
Length	45"	44"
Width	22"	22"
Height (ground to top of seat back)	37"-39"	37"-39"
Height (ground to top of headrest)	43.5"-48.5"	43.5"-48.5"
Ground to Top of Seat	21.75"-23.75"	21.75"-23.75"
Ground to Top of Deck	6.75"	6.75"
Deck to Top of Seat	15.75"-17.75"	15.75"-17.75"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	34.7"	33.8"

Weight of Unit: (assembled)		
Front Half	168 lbs.	174 lbs.
Rear Half	50 lbs.	57 lbs.
Seat w/Arms	36 lbs.	36 lbs.
Batteries	32 lbs.	32 lbs.
	24 lbs./pack	24 lbs./pack
	x 2 = 48 lbs.	x 2 = 48 lbs.

Tires:		
Standard Tires	Solid, low profile	Solid, low profile
Optional Tires	N/A	N/A
Front	9"	9"
Rear	9"	9"
Rear Anti-Tip Wheels	2"	2"

Standard Seat		
Seat Back Height - no headrest	High Back	High Back
Seat Back Height - with headrest	17"	17"
Seat Back Height - with headrest	23"-26"	23"-26"
Width x Depth	17" x 16"	17" x 16"
Color	Black Vinyl	Black Vinyl

Optional Seat:		
Back Height - no headrest	Stadium Style	Stadium Style
Back Height - no headrest	16"	16"
Width x Depth	20" x 18"	20" x 18"
Color	Black Vinyl	Black Vinyl

Warranty:		
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics ²	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

¹ Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain, and battery type.
² Electronics warranty excludes batteries.


Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



1. An Infinite Adjustable Tiller (no knuckle) at your fingertips!



2. The tiller folds easily and locks in place.



3. Ultrabrite LED headlight and taillight!



4. USB & Battery Charging Ports conveniently located on the tiller under the control panel.



5. A stylish new profile offers ample leg room!



Disassembles into 6 pieces!
3 wheel model shown



Phone: 800-265-9615
www.goldentechnologies.ca

Find us on Facebook
www.facebook.com/goldentechcanada