




Standard Color



Black
Carbon Fiber

**ONLY
33 lbs.***
to carry
but holds up to
300 lbs.!

Cricket GP302



The all-new Cricket Power Wheelchair is stylish, agile, and extremely portable. Crafted with an incredibly strong yet lightweight 100% carbon fiber frame, the Cricket weighs as light as 33 lbs.*, making travel a breeze! The slim sporty design is perfect for navigating indoors and out while sightseeing, shopping, or dining. The roomy padded seat cushion and durable mesh back align with the natural curves of your body providing long-lasting comfort while on the go. The Cricket features a well-styled joystick designed to make everyday tasks easier such as driving up to a table or maneuvering through a doorway. It has an impressive 22" turning radius and offers a 300 lb. weight capacity, maximum speed of 3.7 mph, and an operating range of up to 9.3 miles on airline-approved 12 AH 24V lithium-ion batteries. It also offers convenient offboard charging and folds effortlessly for transport. The Cricket offers excellent performance and portability for your active lifestyle!

Standard Features:

- 100% Carbon Fiber Frame
- 300 lb. Weight Capacity
- Removable Compact Joystick with USB Charger
- 35 Amp Controller
- Convenient Under Seat Storage Bag
- Removable Padded Seat
- Offboard Charging Capability
- Easy One-Step Folding
- Stands on End for Space-saving Storage
- Airline-approved Lithium-ion Battery





FOLDABLE POWER CHAIR

| | |
|--|----------------------|
| Model Number | GP302 |
| Class I Medical Device | Yes |
| Weight Capacity | 300 lbs. |
| Drive Wheels | Rear |
| Maximum Speed | 3.7 mph |
| Operating Range ¹ | 9.3 miles |
| Ground Clearance | 2.6" |
| Turning Radius | 22" |
| Type Batteries | 12AH 24V Lithium-Ion |
| Color | Black Carbon Fiber |
| Freewheel Mode | Yes |
| Electronic Speed Control | Yes |
| Electro-Mechanical Brakes | Yes |
| Charger | Off Board |
| Controller | 35 Amp |
| Dimensions: | |
| Length | 36.6" |
| Width | 22.25" |
| Height (ground to top of back of seat) | 36.5" |
| Ground to Top of Seat | 19.3" |
| Number of Seat Height Positions | N/A |
| Size of increments | N/A |
| Weight of Unit: (assembled) | 39.3 lbs. |
| Weight of Unit: (minus parts)* | 33.2 lbs. |
| Front Half | N/A |
| Rear Half | N/A |
| Seat w/Arms | N/A |
| Batteries | 3.8 lbs. |
| Tires: | |
| Drive Wheels (Rear) | 8.5" |
| Casters (Front) | 7" |
| Rear Anti-Tip Wheels | 1.6" |
| Standard Seat: | |
| Back Height - no headrest | 19.3" |
| Width x Depth | 17.3" x 16.5" |
| Color | Black |
| Footrest: | |
| Height adjustable | N/A |
| Angle adjustable | N/A |
| Warranty: | |
| Frame | 5 Years |
| Drive Train | 13 Months |
| Electronics ² | 13 Months |
| Batteries | 12 Months |

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.

* Excludes battery, joystick, cushion and bag.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



**ONLY
33 lbs.***
to carry
but holds up to
300 lbs.!


Sporty, strong and high-performing with a 100% carbon fiber frame!



Compact Joystick



Convenient USB Charging Port



Flip Up Arm Rests

**FOLD
& STOW!**

Compact Dimensions
30.2" x 22.25"
x 14"



East One-Step Folding for Transport



Stands on End for Storage



Removable 1.25" Thick Padded Seat



Phone: 800-265-9615
www.goldentechnologies.ca



www.facebook.com/goldentechcanada

SHT-GP302-CAN-0326