



Standard Colors:



Red



Blue

Compass™ Sport GP605



The Golden Compass™ Sport features a stylish, sleek design built for maneuverability both indoors and out with durable, solid color injected plastic shrouds in red or blue for a sporty new look. It offers a range of comfortable seat sizes available in 16", 18" and 20" standard widths and an upgraded 22" wide seat. All four seat widths come with a luxury high-back Captain's style seat with a headrest for ultimate comfort when on the go. The Dynamic LiNX controller offers a better rider experience using patented Dynamic Load Compensation Technology, which improves the terrain climbing ability of the Compass Sport's springless articulating chassis. Learning to drive the Compass Sport is easier and a more natural experience than before!

Standard Features:

- Solid color injection molded plastics
- Nylon drive wheels for reduced weight and modern look
- Oil-embedded bearings for increased durability
- Silver powder coated arm beams
- Four reflective stickers
- Larger footplate
- Optional 50 Ah batteries for up to 34-mile operating range



Golden Compass™ Sport Specifications

Model	GP605
Model Name	Compass Sport
Class I Medical Device	Yes
Maximum Incline	6°
Weight Capacity	300 lbs
Drive Wheels	Center
Maximum Speed	4 mph
Operating Range ¹	24 miles / up to 34 miles
Ground Clearance	3.25"
Turning Radius	24"
Type Batteries	2-U1 / up to 50 Ah
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp
Dimensions:	
Length	39.5"
Width (outside of tires)	24"
Height (ground to top of back of seat)	40" - 42"
Height (ground to top of back of headrest)	47" - 52.5"
Ground to Top of Seat	19" - 21"
Number of Seat Height Positions	3
Size of Increments	1"
Weight of Unit: (assembled)	179 lbs
Base	97 lbs
Seat with Arms	40 lbs
Batteries	22 lbs each / 30 lbs each
Tires:	
Center	10"
Casters	6"
Standard Seat:	
Back Height - no headrest	23"
Back Height - with headrest	30.5" - 34.5"
Width x Depth	18" x 18" / 20"
Color	Black Vinyl
Optional Seats:	
Back Height - no headrest	19"
Back Height - with headrest	26.5" - 30"
Width x Depth	16", 18", 20", or 22" x 16", 18" or 20"
Color	Black Vinyl
Footrest:	
Height adjustable	Yes
Angle adjustable	2.5" - 4.5"
Angle adjustable	0 - 10 degrees
Warranty:	
Frame	Lifetime
Drive Train	13 Months
Electronics ²	13 Months
Batteries	12 Months

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain & battery type.

² Electronics warranty excludes batteries.



Solid Color Injection Molding Plastics in 2 Gorgeous Colors!



Shown in Blue



Phone: 800-265-9615
www.goldentechnologies.ca

Find us on
Facebook

www.facebook.com/goldentechcanada