



Standard Colors:



Red



Stone Grey



Desert Tan

Eagle 4-Wheel GR595



If you're looking for the ultimate outdoor driving experience, look no further than the new Golden Eagle™. In a class all its own, the Eagle's stylish yet rugged, off-road design is built for high performance and maneuverability with front and rear suspension for smoothly traversing outdoor terrain and a powerful 140 Amp controller for outstanding climbing ability. It offers a full LED lighting package, an adjustable tiller and a luxury high-back Captain's seat for maximum comfort when on the go, plus durable vibrant red color shrouds for a sporty look. With a 500 lb. weight capacity, 13" pneumatic tires, up to 32-mile operating range on optional 75 AH batteries, plus a water-resistant storage compartment that comes standard, the Golden Eagle offers everything you need to enjoy the great outdoors!

Standard Features:

- 500 lb. Weight Capacity
- Integrated LCD Console Display
- Full LED Lighting Package
- Built-in Speaker with Bluetooth® & USB Connectivity
- Lockable Water-resistant Storage Compartment
- Rearview Mirrors
- Cell Phone Holder & USB Charger
- Convenient Drink Holders
- Infinite Position Tiller
- High-back, Limited Recline Captain's Seat
- Large Front Basket
- 13" Pneumatic Tires
- 2-NF22 Batteries

**Upgrade with
75 AH Batteries for
Extended Range!**



Item	GR595
Class I Medical Device	Yes
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	8 mph
Operating Range ¹	23 miles / 32 miles
Ground Clearance	5"
Turning Radius	69"
Type Batteries	2-NF22 / 2-75 AH
Color	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Mechanical Hand Brake	No
Charger	Off Board
Controller	140 Amp
Dimensions:	
Length	56.5"
Width	27"
Height (ground to top of back of seat)	49" - 50"
Height (ground to top of back of headrest)	57" - 62.5"
Ground to Top of Seat	27" - 28"
Ground to Top of Deck	7.25"
Deck to Top of Seat	19.75" - 20.75"
Number of Seat Height Positions	2
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight of Unit: (assembled)	
Base	300 lbs
Seat w/Arms	239 lbs
Batteries	49.5 lbs
	37 lbs ea (NF22)
	52 lbs ea (75 AH)
Tires:	
Front ³	Pneumatic
Rear ³	13"
Standard Seat:	
Slide Seat	Captain's High-Back
Back Height - no headrest	Standard
Back Height with headrest	23.75"
Width x Depth	30.75" - 36.25"
Color	20" x 18"
Warranty:	
Frame	Black
Drive Train	5 Years
Electronics ²	2 Year
Batteries	2 Year
	1 Year

¹ Battery range at 200 lbs., and will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.

³ Inflate tires to 35 psi.



Rear view mirrors, cell phone holder and built-in speaker system with Bluetooth® & USB connectivity.



Integrated LCD console display monitors speed, odometer, fault codes & battery life.



Lockable water-resistant storage compartment comes standard and attaches to the accessory bracket.



Ultrabright LED rear lights.



Phone: 800-265-9615
www.goldentechnologies.ca

Find us on Facebook
www.facebook.com/goldentechcanada