



Standard Colors:



Crimson Red



Stone Grey

Standard Features:

- 500 lb. Weight Capacity
- LCD Color Console Screen with 3 Customizable Gauge Layouts
- Turning Speed Reduction Sensor
- Front & Rear Adjustable Suspension
- Powerful 200 Amp Controller
- Full LED Lighting Package
- Bluetooth® Speakers
- Lockable Water-resistant Storage Compartment
- Rear View Mirrors
- Cell Phone Holder
- USB Charging Port
- Convenient Drink Holders
- Infinite Position Tiller
- High-back, Limited Recline Captain's Seat
- Width, Angle & Height Adjustable Armrests
- Large Front Basket
- 12" Pneumatic Tires
- NFC Key Tag & Conventional Key Ignition
- 2-75 AH Batteries

Eagle A/T 4-Wheel GR596



Take your outdoor driving experience to new heights with the Golden Eagle All-Terrain™ scooter. In a class all its own, the Eagle A/T's stylish yet rugged, off-road design is built for superior performance and maneuverability with front and rear adjustable suspension for smoothly traversing outdoor terrain and a powerful 200 Amp controller for outstanding climbing ability. It offers full LED lighting with brake light and turn signals, an infinitely adjustable wraparound tiller, and a reclining, high-back Captain's seat for maximum comfort, plus vibrant red or grey color shrouds for a sporty look. With a 500 lb. weight capacity, 12" pneumatic turf tread tires, up to 34-mile operating range on 75 AH batteries, plus a water-resistant rear carry case for personal items, the Golden Eagle A/T offers everything you need to enjoy your next adventure!



Upgrade to
100 AH Batteries for
Extended Range!

Item	GR596
Class I Medical Device	Yes
Weight Capacity	500 lbs
Drive Wheels	Rear
Maximum Speed	11 mph
Operating Range ¹	34 miles w/ 75 AH Batteries
Ground Clearance	4"
Turning Radius	76"
Type Batteries	75 AH - 100 AH
Color	Red & Grey
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Mechanical Hand Brake	No
Charger	Off Board
Controller	200 Amp
Dimensions:	
Length	63"
Width	31.7"
Height (ground to top of back of seat)	56.5"
Ground to Top of Seat	26.5"
Ground to Top of Deck	7"
Deck to Top of Seat	19.5"
Number of Seat Height Positions	2
Size of Increments	1"
Front Axle to Rear Axle	40"
Weight of Unit: (assembled)	377 lbs w/ 75 AH 407 lbs w/ 100 AH
Base	208 lbs
Seat w/Arms	54 lbs
Batteries (each)	52 lbs (75 AH) 67 lbs (100 AH)
Tires:	
	Pneumatic
Front ³	13" x 6.5"
Rear ³	13" x 6.5"
Standard Seat:	
Slide Seat	Captain's High-Back
Seat Back Height	Standard
Width x Depth	23.75"
Color	20" x 18"
Warranty:	
Frame	Black Vinyl
Drive Train	5 Years
Electronics ²	2 Years
Batteries	2 Years
	1 Year

¹ Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.

³ Inflate tires to 28 psi.



High-Tech NFC Key Tag & Conventional Key Ignition



LCD Color Console Screen with 3 Customizable Gauge Layouts



Lockable water-resistant storage compartment comes standard and attaches to the accessory bracket.



Ultrabright LED rear lights.



Phone: 800-265-9615
www.goldentechnologies.ca

Find us on
Facebook

www.facebook.com/goldentechcanada