




## Standard Fabrics:



Antler



Lakefront



NEW! Umber

The Imagine contemporary plush fabric pattern comes in beautiful, rich colors. It features our exclusive SofTouch™ with Stain Defense coating.

The Sinda polyurethane pattern is soft and supple and was designed to be breathable yet highly durable. It provides a realistic leather-like look with a subtle texture. The composite material and polyurethane resin coating make this fabric water resistant, flexible, and lightweight.

Crossover fabrics are unavailable.

**Elara**  
PL/PR118



*Comfort Zone*™  
by GOLDEN

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Elara is offered in two sizes with a back style that provides maximum comfort in a tasteful seam design. It features a headrest, supportive chaise pad, and a USB port on the hand control. Available in luxurious new fabrics to match any décor. Enjoy the powerful energy and take time to better connect with yourself in your DeLuna recliner by Golden!

## Standard Features:

- Easy Grip Armrest
- Soft-Feel Hand Control
- USB Charging Port
- SmartTek® Diagnostic System
- New! Sinda Umber Urethane Fabric. Its water-resistant and wipeable surface ensures easy cleaning and superior wear resistance, all backed by a 5-year warranty



# DeLUNA Series PL/PR118

Model Number	118-MSM	118-LAR
Size	Med-Small	Large
Class II Medical Device	YES	YES
Weight Capacity	400 LBS	400 LBS
Positioning	3	3
Comfort Zone(s)	1	1
Twilight Technology	NO	NO
ZG+	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	YES
Back Type	Seam	Seam
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	Standard	Standard
Arm Rest Covers	NO	NO
Overall Width	32"	34"
Overall Height	43"	46"
Seat Width	21"	23"
Floor to Top of Seat	18.5"	19.5"
Seat Depth	21"	21"
Seat to Top of Back	29"	31"
Seat to Top of Arm	7"	6"
Distance Required from Wall Reclined	19"	20"
Extension Reclined from Back of Seat to End of Foot Rest	42"	42"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Arm	Inside Arm
Extra Pocket Option	NO	NO
Heat & Vibrating Massage Option	NO	NO
HeatWave Technology Option	NO	NO
Articulating Headrest Option	NO	NO
Power Adjustable Lumbar Option	NO	NO
Foot Rest Extension Option	NO	NO
Straight Lift Option	Standard	Standard
USB Charging Port	Hand Control	Hand Control
Battery Backup	9 volt	9 volt
SmartTek Diagnostic System	Standard	Standard
Place of Origin	Imported	Imported



Elegant seamed back design and supple urethane fabric make this chair both pleasing to the eye and highly durable.

Shown in Lakefront



Shown in Lakefront

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.

## Comfort Zone™ by GOLDEN

**Comfort Zone 1** provides ergonomic placement of your lower body allowing your muscles to relax for long-lasting comfort. Comfort Zone 1 can be adjusted using the two buttons on the hand control.



LOWER BODY



UPPER BODY



HEADREST



LUMBAR



TWILIGHT TILT

**WE PATENTED COMFORT™**



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-265-9615  
www.goldentechnologies.ca

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with  
UL 73 and  
CSA C22.2 No. 68



www.facebook.com/goldentechcanada

SHT-PR118-CAN-0326