





Standard Fabrics:




Calypso



Camel



Hazelnut



Shiraz



Sterling

This fabric has a modern, updated look with rich color options. It features our exclusive SofTouch™ with Stain Defense coating.

* An upcharge will apply for all crossover fabrics.

Standard Features:

- Patented Twilight Tilt Technology
- Patented MaxiComfort® Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable AutoDrive™ Hand Control
- USB Charging Port
- Extra Pocket
- SmartTek® Diagnostic System

Upgrade Options:

- Deluxe Therapeutic Heat & Vibrating Massage

Cloud
PR515



Comfort Zone™
by GOLDEN

The Cloud™ is a luxury power lift recliner offering our unique, ergonomic seating system for long-term comfort and support. It features our custom bucket seat and plush backrest. Positioning options are endless with our 5-zone positioning system, which includes our patented Twilight Technology, patented MaxiComfort® Positioning, adjustable headrest and lumbar support. Healthy benefits include being able to elevate the feet above the heart while experiencing a feeling of weightlessness. Now available in two sizes!



MAXICOMFORT® Series PR515

Model Number	PR515-SME	PR515-MLA
Size	Small/Med	Med/Large
Class II Medical Device	YES	YES
Weight Capacity	400 LBS	400 LBS
Positioning	MaxiComfort	MaxiComfort
Comfort Zone(s)	5	5
Twilight Technology	YES	YES
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	NO
Back Type	Biscuit	Biscuit
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	YES	YES
Arm Rest Covers	NO	NO
Overall Width	35.5"	39.5"
Overall Height	44"	46.5"
Width Between Arms	20"	20"
Floor to Top of Seat	20"	21.5"
Seat Depth	18.5"	20"
Seat to Top of Back	28"	31"
Seat to Top of Arm	8"	8.5"
Distance Required from Wall Reclined	23"	24"
Extension Reclined from Back of Seat to End of Foot Rest	43"	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Pocket	Inside Pocket
Extra Pocket Option	Standard	Standard
Heat & Massage Option	YES	YES
HeatWave Technology Option	YES	YES
Nirvana with HeatWave Option	NO	NO
Articulating Headrest Option	Standard	Standard
Adjustable Lumbar Option	Standard	Standard
Foot Rest Extension Option	YES	YES
Straight Lift Option	NO	NO
USB Charging Port	Control Box	Control Box
Battery Backup	YES	YES
SmartTek Diagnostic System	Standard	Standard
Place of Origin	USA	USA



TWILIGHT

The World's Most Advanced
Power Positioning Lift Recliner
with Patented Twilight Technology!



Shown in
Coffee Bean

Healthy Benefits:

- Elevates Feet Above Heart
- Alleviates Pressure On Lungs
- Distributes Weight
- Relieves Pressure Points
- Encourages Proper Posture
- Supports Natural 'S' Shape of Spine
- Available in Luxurious Brisa® Fabrics With Breathable Technology

Comfort Zone™

by GOLDEN

Comfort Zones 1 - 5 This positioning system incorporates all five adjustable comfort zones to create the ultimate wellness comfort experience with our patented Twilight Technology. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive™ hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-265-9615
www.goldentechnologies.ca



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

Find us on
Facebook

www.facebook.com/goldentechcanada

LIT-SHT-PR515 CAN-0326