



Austin UC342



Comfort Zone[™]
by GOLDEN

Take a vacation from your busy day with the Austin by Golden. This power lift recliner offers a plush seamed back and supportive chaise pad providing head to toe comfort you can really sink into. And for those days when tired muscles need some extra TLC, simply turn on the heat and vibrating massage for an enhanced relaxation experience!

Standard features:

- Heat & Vibrating Massage
- SmartTek Diagnostic System

Standard Fabrics:



Graphite



Maple



UC342 - Austin

Model Number	UC342
Size	Med
Weight Capacity	400 LB
Positioning	3
Comfort Zone(s)	1
Eclipse Technology	No
Chaise Seat	Yes
Pocketed Coil Seat	No
Back Type	Seam
Removable Back Type	Bracket
Overall Width	32"
Overall Height	41"
Overall Depth (Seated)	35.5"
Width Between Arms	21"
Floor to Top of Seat	19"
Seat Depth	21"
Seat to Top of Back	28"
Seat to Top of Arm	6"
Distance Required from Wall Reclined	17"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Lifetime Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	No
Warranty on Mechanical Labor	1 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric and Foam	1 YR
Head Rest Cover	No
Arm Covers	No
Straight Lift	Yes
Hand Control Routing	Inside Arm
Extra Pocket Available	Standard
Heat & Massage Available	2-Motor Standard
HeatWave Technology Available	No
Articulating Headrest Available	No
Adjustable Lumbar Available	No
Foot Rest Extension Available	No
USB Charging Port	No
Battery Backup	Yes
Place of Origin	China



Smart Tek DiagnosticSystem STANDARD



Comfort Zone[™]
by GOLDEN



WE PATENTED COMFORT[™]

 **GOLDEN**

Designed and engineered in Old Forge, PA. Assembled in China.

Literature is current at time of printing. Golden reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

Telephone : 800-265-9615
www.goldentechnologies.ca

Like us on Facebook
Facebook.com/GoldentechCanada

SHT-UC342-CAN-0326