



Standard Colors:



Red



Blue

BuzzAbout GP164

BUZZABOUT

Enjoy full-time luxury on the go in the new BuzzAbout! Perfect for travel, the BuzzAbout power wheelchair is in a class all its own in our popular BuzzAround series of mobility products. The BuzzAbout has a stylish captain's seat with contoured support and higher back for longer seating tolerance and comfort. It features motion sensing LED running lights on the footplate and rear handle activated when the joystick is moved. It also offers convenient offboard charging and easy disassembly for transport. The BuzzAbout offers excellent performance and maneuverability indoors and out for your active lifestyle.

Features:

- Motion Sensing LED Running Lights on Footplate & Rear Handle
- Quick & Easy Disassembly into Four Pieces
- Offboard Charging Capability
- 300 lb. Weight Capacity



Item

Description	GP164
Class I Medical Device	Yes
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range up to ¹	11 miles
Under Deck Clearance	2.5"
Turning Radius	28.5"
Type Batteries	2-22AH
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off-board
Controller	40 Amp Remote

Dimensions

Length	36"
Width	23" at front casters
Height (ground to top of back of seat)	36.75-38.75"
Ground to Top of Seat	21.5"-23.5"
Number of Seat Height Positions	3
Size of Increments	1"

Weight of Unit: (assembled)

Front Half	29 lbs.
Rear Half	33 lbs.
Seat w/Arms	34 lbs.
Batteries	28 lbs.

Tires:

Front	N/A
Rear	9"
Casters	6"
Rear Anti-Tip Wheels	2"

Standard Seat

Slide Seat	N/A
Back Height (no headrest)	18"
Width x Depth	17" x 16"
Color	Black Vinyl

Footrest:

Height adjustable	No
Angle adjustable	Yes

Warranty:

Frame	5 Years
Drive Train	13 Month Warranty
Electronics ²	13 Month Warranty
Batteries	1 Year Warranty

¹. Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

². Electronics warranty excludes batteries.



Front Motion Sensing LED Light



Rear Motion Sensing LED Light



GP164 Shown Disassembled



Phone: 800-265-9615
www.goldentechnologies.ca



www.facebook.com/goldentechcanada