




## Standard Color



Black

AS  
LIGHT AS  
**52.6**  
lbs.!\*

## Stride GP301



The all-new Stride Power Wheelchair will take you where you want to go in comfort and style! The Stride is stylish with a two-tone grey and black seat with removable cushion and a tall back for long-lasting comfort while on the go. Modern-looking front and back comfort spring suspension add to the smooth ride. The sporty design is perfect for navigating indoors and out while sightseeing, shopping, or dining. This portable folding power chair is lightweight with an aluminum frame and stows easily in the trunk of your car. Weighing as light as 52.6 lbs.\*, the Stride is the perfect choice to take on road trips, cruise ships and airplanes. The Stride features a well-styled joystick designed to make everyday tasks easier such as driving up to a table or maneuvering through a doorway. It has an impressive 23.6" turning radius and offers a 265 lb. weight capacity, maximum speed of 3.7 mph, and an operating range of up to 9.3 miles on airline-approved 10AH 24V lithium-ion batteries. It also offers convenient offboard charging and folds effortlessly for transport.

## Standard Features:

- Aluminum Frame
- 265 lb. Weight Capacity
- Removable Compact Joystick with USB Charger
- 50 Amp Controller
- Comfort Spring Suspension
- Removable Padded Seat
- Convenient Under-Seat Storage & Back Seat Pocket
- Offboard Charging Capability
- Airline-approved Lithium-ion Battery
- Easy One-Step Folding
- Stands on End for Storage





## FOLDABLE POWER CHAIR

|                                        |                                        |
|----------------------------------------|----------------------------------------|
| <b>Model Number</b>                    | <b>GP301</b>                           |
| Class I Medical Device                 | Yes                                    |
| Weight Capacity                        | 300 lbs.                               |
| Drive Wheels                           | Rear                                   |
| Maximum Speed                          | 3.7 mph                                |
| Operating Range <sup>1</sup>           | 9.3 miles                              |
| Ground Clearance                       | 1.6"                                   |
| Turning Radius                         | 23.6"                                  |
| Type Batteries                         | 10AH 24V Lithium-Ion                   |
| Color                                  | Black                                  |
| Freewheel Mode                         | Yes                                    |
| Electronic Speed Control               | Yes                                    |
| Electro-Mechanical Brakes              | Yes                                    |
| Charger                                | 2 Amp, Off Board                       |
| Controller                             | 50 Amp                                 |
| <b>Dimensions:</b>                     |                                        |
| Length                                 | 42.5"                                  |
| Width                                  | 24.4"                                  |
| Height (ground to top of back of seat) | 38.6"                                  |
| Ground to Top of Seat                  | 22"                                    |
| Number of Seat Height Positions        | N/A                                    |
| Size of Increments                     | N/A                                    |
| <b>Weight of Unit: (assembled)</b>     | <b>58.0 lbs.</b>                       |
| <b>Weight of Unit: (minus parts)*</b>  | <b>52.6 lbs.</b>                       |
| Front Half                             | N/A                                    |
| Rear Half                              | N/A                                    |
| Seat w/Arms                            | N/A                                    |
| Batteries                              | 3.8 lbs                                |
| <b>Tires:</b>                          |                                        |
| Drive Wheels (Rear)                    | 12"                                    |
| Casters (Front)                        | 8"                                     |
| Rear Anti-Tip Wheels                   | 2"                                     |
| <b>Standard Seat:</b>                  |                                        |
|                                        | Sling Back with Removable Seat Cushion |
| Back Height - no headrest              | 16"                                    |
| Width x Depth                          | 17" x 18"                              |
| Color                                  | Black                                  |
| <b>Footrest:</b>                       |                                        |
| Height adjustable                      | N/A                                    |
| Angle adjustable                       | N/A                                    |
| <b>Warranty:</b>                       |                                        |
| Frame                                  | 5 Years                                |
| Drive Train                            | 13 Months                              |
| Electronics <sup>2</sup>               | 13 Months                              |
| Batteries                              | 12 Months                              |

<sup>1</sup> Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

<sup>2</sup> Electronics warranty excludes batteries.

\* Excludes battery, joystick, cushion and bag.

**AS  
LIGHT AS  
52.6  
lbs.!\***



Portable & high-performing  
with an aluminum frame!



East One-Step Folding with Velcro  
Closure Strap for Transport

**FOLD  
& STOW!**  
Compact Dimensions  
42.5" x 24.4"  
x 15.7"

Stands On End  
for Storage



Comfort Spring Suspension for a  
Smooth Ride



Removable Padded Seat



Phone: 800-265-9615  
[www.goldentechnologies.ca](http://www.goldentechnologies.ca)

**f** Find us on  
**Facebook**

[www.facebook.com/goldentechcanada](http://www.facebook.com/goldentechcanada)

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

SHT-GP301-CAN-0326