



Cloud PR527




Comfort Zone™
by GOLDEN

The Cloud™ luxury power lift recliner blends classic design with innovative features like an ultra-plush, overstuffed bucket-back and a uniquely contoured memory foam bucket seat that cradles your hips and encourages healthy posture. Positioning options are endless with our 5-zone positioning system, which includes our patented Twilight Technology, patented MaxiComfort® positioning, adjustable headrest, and lumbar support. Upgrade your comfort and wellness experience by adding our Nirvana™ Swedish Air Massage with full-body, Far Infrared HeatWave™ Technology or Deluxe Heat and Vibrating Massage.




FDA Class I Medical Devices designed by Golden Technologies are intended to assist persons with chronic physical disabilities.

Standard Fabrics:



Ink




Sand



Slate



Spa



Stone

This fabric has a modern, updated look with rich color options.

* An upcharge will apply for all crossover fabrics.

Standard Features:

- Patented Twilight Tilt Technology
- Patented MaxiComfort® Positioning with Four-Way Motion Seat
- Adjustable Headrest
- Adjustable Lumbar
- Programmable AutoDrive™ Hand Control
- USB Charging Port
- Extra Pocket
- SmartTek® Diagnostic

Optional Features*:

- HeatWave™ Technology
- Deluxe Heat & Vibrating Massage
- Nirvana™ Swedish Air Massage with Far Infrared Heat Therapy

*Additional lead time may apply

MAXICOMFORT® Series PR527

Model Number	PR527-SME	PR527-MLA
Size	Small/Med	Med/Large
Class I Medical Device	YES	YES
Weight Capacity	400 LBS	400 LBS
Positioning	Twilight	Twilight
Comfort Zone(s)	5	5
Twilight Technology	YES	YES
ZG+ Position	YES	YES
Chaise Seat	YES	YES
Pocketed Coil Seat	NO	NO
Back Type	Seam	Seam
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	NO	NO
Arm Rest Covers	NO	NO
Overall Width	38.5"	38.5"
Overall Height	44"	47"
Seat Width	23.5"	23.5"
Floor to Top of Seat	19.5"	22.5"
Seat Depth	20.5"	32"
Seat to Top of Back	29"	28"
Seat to Top of Arm	8.5"	8.5"
Distance Required from Wall Reclined	21"	21.5"
Extension Reclined from Back of Seat to End of Foot Rest	42"	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Outside Arm	Outside Arm
Extra Pocket Option	Standard	Standard
Heat & Vibrating Massage Option	YES	YES
HeatWave Technology Option	YES	YES
Nirvana with HeatWave Option	YES	YES
Articulating Headrest Option	Standard	Standard
Power Adjustable Lumbar Option	Standard	Standard
Foot Rest Extension Option	NO	NO
Straight Lift Option	Standard	Standard
USB Charging Port	Right Arm	Right Arm
Battery Backup	Lithium-ion	Lithium-ion
SmartTek Diagnostic System	Standard	Standard
Place of Origin	USA	USA

TWILIGHT

The World's Most Advanced
Power Positioning Lift Recliner
with Patented Twilight Technology!



Shown in Slate

Healthy Benefits:

- Elevates Feet Above Heart
- Supports Natural 'S' Shape of Spine
- Alleviates Pressure On Lungs
- Available in Luxurious Brisa® Fabrics With Breathable Technology
- Distributes Weight
- Relieves Pressure Points
- Encourages Proper Posture




Comfort Zone™ by GOLDEN

Comfort Zones 1 - 5 This positioning system incorporates all five adjustable comfort zones to create the ultimate wellness comfort experience with our patented Twilight Technology. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive™ hand control.



WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-265-9615
www.goldentechnologies.ca



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

Find us on
Facebook

www.facebook.com/goldentechnologies

LIT-SHT-PR527-CAN0326