





Standard Fabrics:



Bozeman Cub



Bozeman Mineral

This luxurious wipeable fabric has a modern, updated look with rich color options.

Crossover fabrics are unavailable.

Aspen
PL104



Comfort Zone™
by GOLDEN

Inspired by the moon, the DeLuna™ Series of luxury power lift recliners is designed to offer you the most calming comfort and restorative sleep. The Aspen has a modern look with a plush seams back and supportive chaise pad providing head-to-toe comfort. Stylish front-load open arms add modern functionality and accessibility, featuring an integrated wireless smartphone charger, conveniently located USB and USB-C charging ports, and a built-in drink holder, keeping everything you need comfortably within reach. It features a power adjustable headrest and footrest extension for lasting comfort. The Aspen is offered in a luxurious wipeable fabric in colors to match any décor. Enjoy the powerful energy and take time to better connect with yourself in your DeLuna recliner by Golden!



Standard Features:

- Power Adjustable Headrest
- Footrest Extension
- Wireless Smartphone Charger
- USB & USB-C Charging Ports
- Built-in Drink Holder
- Extra Pocket
- SmartTek® Diagnostic System



DeLUNA Series PL104



Model Number	PL104-MED	PL104-LAR
Size	Medium	Large
Weight Capacity	400 LBS	400 LBS
Positioning	3	3
Comfort Zone(s)	2	2
Twilight Technology	NO	NO
Chaise Seat	YES	YES
Pocketed Coil Seat	YES	YES
Back Type	Seam	Seam
Removable Back Type	Quick Release	Quick Release
Head Rest Cover	NO	NO
Arm Rest Covers	NO	NO
Overall Width	39"	41"
Overall Height	43"	47.5"
Seat Width	21"	23"
Floor to Top of Seat	19.5"	20"
Seat Depth	21"	23"
Seat to Top of Back	29"	33"
Seat to Top of Arm	8"	8"
Distance Required from Wall Reclined	15"	17.5"
Extension Reclined from Back of Seat to End of Foot Rest	45.5"	46.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES
Warranty on Electrical Parts	3 YR	3 YR
Pro Rated Years 4 Thru 7	YES	YES
Warranty on Mechanical Labor	3 YR	3 YR
Warranty on Electrical Labor	1 YR	1 YR
Hand Control Routing	Inside Pocket	Inside Pocket
Extra Pocket Option	Standard	Standard
Heat & Vibrating Massage Option	NO	NO
HeatWave Technology Option	NO	NO
Nirvana with HeatWave Option	NO	NO
Adjustable Headrest Option	Standard	Standard
Adjustable Lumbar Option	NO	NO
Foot Rest Extension Option	Standard	Standard
Straight Lift Option	NO	NO
USB Charging Port	Right Arm	Right Arm
Battery Backup	YES	YES
Place of Origin	Thailand	Thailand




Shown in Bozeman Cub



Power Headrest



Footrest Extension



Wireless Phone Charger



Drink Holder

Designed and engineered in Old Forge, PA. Assembled in Thailand using Golden components.



Comfort Zone™
by GOLDEN

THE HALO provides ergonomic placement of your lower body allowing muscles to relax, while a power adjustable headrest and footrest extension refine long-lasting comfort with minimal effort.



WE PATENTED COMFORT™



Phone: 800-265-9615
www.goldentechnologies.ca

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

Find us on
Facebook

www.facebook.com/goldentechcanada

LIT-SHT-PL104-CAN-0926