



Standard Colors:



Red



Blue

Buzzaround 3 wheel GB121B-SHZ 4 wheel GB124A-SHZ



The next generation of the popular Buzzaround XLSDHD features compact portability with a heavy-duty weight capacity! The XLSDHD features a robust 325 lb. weight capacity and our Comfort Spring Suspension for a smoother ride. It offers improved standard features like a larger wraparound Delta tiller, higher ground clearance of 4", and an operating range of up to 9 miles on the standard 12 AH batteries. Upgrade with the new XR Extended Range Battery option for even more mileage! The Buzzaround XLSDHD also offers a new stylish, contoured stadium style seat for more comfort and increased seating tolerance. This scooter disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!

Features:

- 325 lb. Weight Capacity
- Comfort Spring Suspension
- Manually Adjustable Preload on Rear Suspension
- Improved Drivability
- Higher Ground Clearance
- LED Angle-adjustable Headlight
- USB Port on the Tiller
- Phone & Cup Holders
- Quick & Easy Disassembly

Upgrade with NEW XR Extended Range Battery Option!

(Packaged separately as an accessory but ships with the unit.)





LIGHTWEIGHT HEAVY-DUTY PORTABLE SCOOTER

Item	GB121B-SHZ	GB124A-SHZ
Description	XLSHD	XLSHD
Medicare Code	Not Coded	Not Coded
Weight Capacity	325 lbs.	325 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to ¹	9 miles/18 miles	9 miles/18 miles
Under Deck Clearance	4"	4"
Turning Radius	37.5"	44"
Type Batteries	2-12 AH / 2-22 AH	2-12 AH / 2-22 AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp	50 Amp

Dimensions

Length	41"	40"
Width	20"	20"
Height (ground to top of back of seat)	34.25"-35.75"	34.25"-35.75"
Ground to Top of Seat	22.25"-23.75"	22.25"-23.75"
Ground to Top of Deck	5.5"	5.5"
Deck to Top of Seat	16.5"-18"	16.5"-18"
Number of Seat Height Positions	3	3
Size of Increments	.75"	.75"
Front Axle to Rear Axle	31.5"	30.5"

Weight of Unit: (assembled)³

Front Half	100 lbs.	105 lbs.
Rear Half	35 lbs.	40 lbs.
Seat w/Arms	28 lbs.	28 lbs.
Batteries	21 lbs.	21 lbs.
Batteries	16 lbs. / 30 lbs.	16 lbs. / 30 lbs.

Tires:

	Solid, low profile	Solid, low profile
Front	7.5"	7.5"
Rear	8"	8"
Rear Anti-Tip Wheels	2"	2"

Standard Seat

	Stadium Style	Stadium Style
Seat Back Height - no headrest	14.75"	14.75"
Width x Depth	17" x 14.5"	17" x 14.5"
Color	Black Vinyl	Black Vinyl

Warranty:

Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics ²	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

¹ Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain, and battery type.

² Electronics warranty excludes batteries.



GB124A-SHZ 4 Wheel



Disassembles into 5 Pieces!

3 Wheel Model Shown



Comfort Spring Suspension

Manually adjustable preload on rear suspension!



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com



www.facebook.com/goldentechnologies