



## Buzzaround GP130



The new Buzzaround CarryOn folding power wheelchair is truly in a class all by itself! This unique power chair is super portable with a 21" turning radius and center wheel drive for optimal coverage and maneuverability making it the perfect partner for your next adventure. Our exclusive PowerLink™ Connectivity makes assembly and disassembly a snap! Independent 4-wheel suspension and a springless articulating chassis provide a smooth ride indoors and out. The CarryOn features a 300 lb. weight capacity, front and rear light package, and textured soft touch armrests. Enjoy long range coverage up to 12.5 miles on the airline friendly 12.5 AH lithium-ion battery, which can be removed and taken as a carry on while traveling by air, sea or train. The heaviest piece is 46 lbs. (without seat and batteries) and the height when folded is only 10" for small space storage tucked out of the way. Carry on with life in the CarryOn power chair!



### Standard Colors:



Velvet Red  
(Installed)



Coastal Blue  
(Included in Box)

### Optional Colors:



Crisp White



Canary Yellow



Stone Grey



Orchid Pink



Army Green  
(Matte)



Desert Tan

### Only 46 Lbs.!!\*

\*Without Seat & Battery

**Most Portable, Foldable,  
Lightweight Power Chair with  
Center Wheel Drive!**

**Tightest Turning Radius!**

### Standard Features:

- Center Wheel Drive for Optimal Control & Maneuverability
- Independent 4-Wheel Suspension
- Springless Articulating Chassis
- USB Port
- Front & Rear Lights
- Width Adjustable Armrests
- Textured Soft Touch Armrests
- Angle Adjustable Foot Plate
- Easy Access, Airline-Friendly, 12.5 AH Lithium-ion Battery with Travel Case
- PowerLink™ Connectivity for Easy Disassembly/Assembly
- Removable Storage Basket
- Easily Interchangeable Color Shrouds



## FOLDING POWER WHEELCHAIR

|                                        |                     |
|----------------------------------------|---------------------|
| <b>Model Number</b>                    | <b>GP130</b>        |
| Class I Medical Device                 | Yes                 |
| Weight Capacity                        | 300 lbs.            |
| Drive Wheels                           | Center              |
| Maximum Speed                          | 3.85 mph            |
| Operating Range <sup>1</sup>           | 12.5 miles          |
| Ground Clearance                       | 2"                  |
| Turning Radius                         | 21"                 |
| Type Batteries                         | 12.5 AH Lithium-Ion |
| Colors                                 | Red, Blue           |
| Freewheel Mode                         | Yes                 |
| Electro-Mechanical Brakes              | Yes                 |
| Charger                                | Off Board, 2A       |
| Controller                             | 40 Amp              |
| <b>Dimensions:</b>                     |                     |
| Length                                 | 38"                 |
| Width                                  | 24"                 |
| Height (ground to top of back of seat) | 34.25"              |
| Ground to Top of Seat                  | 21"                 |
| <b>Weight of Unit: (assembled)</b>     | <b>74 lbs.</b>      |
| Base                                   | 46 lbs.             |
| Seat                                   | 22 lbs.             |
| Battery                                | 6 lbs.              |
| <b>Tires:</b>                          |                     |
| Drive Wheels                           | 7.75"               |
| Casters                                | 5.5"                |
| <b>Seat:</b>                           |                     |
| Back Height - no headrest              | 15.5"               |
| Width x Depth                          | 18" x 17"           |
| Color                                  | Black               |
| <b>Warranty:</b>                       |                     |
| Frame                                  | 5 Years             |
| Drive Train                            | 13 Months           |
| Electronics <sup>2</sup>               | 13 Months           |
| Batteries                              | 12 Months           |

<sup>1</sup>. Battery range at 200 lbs., and will vary due to rider weight, drive surface, terrain, and battery type.

<sup>2</sup>. Electronics warranty excludes batteries.



Battery & battery travel bag standard



Folds to only 10" high for storing in small spaces

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Phone: 800-265-9615  
www.goldentechnologies.ca

Find us on  
**Facebook**  
[www.facebook.com/goldentechcanada](http://www.facebook.com/goldentechcanada)

SHT-GP130-CAN-0326