



Standard Colors:



CHERRY RED



ENVY GREEN



PARADISE BLUE



SATIN SILVER

Disassembles Quickly & Easily!

Literider Envoy GP162

LiteRider^{GOLDEN}
ENVY

The LiteRider Envoy Power Wheelchair is both stylish and durable with a tubular design that is unique to the Golden LiteRider family! A Dynamic LiNX controller comes standard and the traditional style makes driving simple and reliable. The new premium seat offers contoured support with an increased back height of 17" for a longer seating tolerance and added comfort. The Envoy, featuring 2-22 amp hour batteries for an operating range of up to 15.5 miles, is Golden's Retail Cash Power Wheelchair! It has a larger footplate, an impressive turning radius of 28.5", and disassembles quickly and easily for transport. When you are looking for a compact power wheel chair that is fun to drive and easy to transport - look no further than the Golden LiteRider Envoy!

Optional Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Pack N' Go
- Seat Belt


Some accessories require a mounting bracket or mounting clips.



Model Number	GP162
Class I Medical Device	Yes
Weight Capacity	300 lbs.
Drive Wheels	Rear
Maximum Speed	3.5 mph
Operating Range ¹	15.5 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Type Batteries	2-22 AH
Colors	Cherry Red, Paradise Blue, Solar Flare Orange, Envy Green, Satin Silver, Sunburst Yellow
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Dynamic Linx
Dimensions:	
Length	36"
Width	22.5"
Height (ground to top of back of seat)	34" - 35"
Height (ground to top of back of headrest)	N/A
Ground to Top of Seat	21.5" - 22.5"
Deck to Top of Seat	N/A
Number of Seat Height Positions	2
Size of Increments	1"
Weight of Unit: (assembled)	123 lbs.
Front Half	37 lbs.
Rear Half	36 lbs.
Seat w/Arms	23 lbs.
Batteries	27 lbs/pack
Tires:	Flat Free
Drive Wheels	9"
Casters	6"
Anti-Tip Wheels	3"
Standard Seat:	Stadium Style
Back Height - no headrest	16"
Back Height with headrest	N/A
Width x Depth	17" x 17"
Color	Black Vinyl
Optional Seats:	Highback Premium
Back Height - no headrest	16" 17"
Back Height with headrest	N/A N/A
Width x Depth	20" x 18" 17" x 16"
Color	Black Vinyl Black Vinyl
Seat Options:	
Pneumatic Rise Seat Option	N/A
Electric Rise Seat Option	N/A
Footrest:	Yes
Height adjustable	Yes
Angle adjustable	Yes
Warranty:	
Frame	5 Years
Drive Train	13 Months
Electronics ²	13 Months
Batteries	12 Months

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.



Watch our video to learn about the benefits of the LiNX controller!

Scan here with smartphone camera.



Disassembles easily for transport!



Ergonomic LiNX controller



This model shown with Premium Seat Option!



Perfect for transport!