




Standard Colors



White



Black Carbon
Fiber Pattern



**Power Folding
Drives Remotely
Tight Turning Radius**

Standard Features:

- Smart Control Mobile App
- Power Folding
- Memory Foam Seating
- Spacious Leg & Foot Area
- Front Suspension
- Omni-Directional Wheels
- Airline Friendly, Lithium-Ion Battery
- System Monitoring Voice Alerts
- Side & Rear Lighting



**Visit App Store
for Remote
Control**

**Exclusive
Omni-Directional
Wheels with
Self-Cleaning
V-Groove!**

**AllySport
GP305**

**GOLDEN
AllySPORT™**

Experience freedom with the new luxury full-time foldable Golden Ally Sport™ power wheelchair. In a class all its own, it features a fashion-forward, sleek design that applies robotics technology for power folding and unfolding and ergonomic seating for maximum comfort. High-tech styling includes built-in side and rear lighting for safety. Designed for seamless navigation on diverse surfaces, the 10" omni-directional front wheels with our exclusive self-cleaning V-groove offer more leg and foot room while gliding effortlessly in any direction with a minimal turning radius. It offers a steady and smooth ride, compliments of its large, pneumatic rear drive wheels and dual front suspension system. Take control through our Smart Control Mobile App, allowing you to manage your device remotely for a personalized and convenient experience. With a long-lasting lithium-ion battery and user-friendly interface, travel anywhere with the Ally Sport – your perfect mobility *ally*!



Model Number	GP305
Class I Medical Device	Yes
Weight Capacity	265 lbs.
Drive Wheels	Rear
Maximum Speed	4 mph
Operating Range ¹	9.5 miles
Ground Clearance	2.5"
Turning Radius	30"
Type Batteries	10 AH Lithium-Ion
Colors	White, Black Carbon Fiber Pattern
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Remote

Dimensions:

Length	42"
Width	25"
Width Between Arms	20" max.
Height (ground to top of back of seat)	37"
Ground to Top of Seat	22.5"

Weight of Unit: (assembled)	92 lbs.
Heaviest Piece	87.5 lbs.
Battery	5 lbs.

Tires:

Drive Wheels	10" Pneumatic
Casters	10" Omni-Directional
Anti-Tip Wheels	2"

Standard Seat:

Back Height - no headrest	18"
Width x Depth	16.5" x 16.5"
Color	Grey or Black

Footrest:

Height adjustable	No
-------------------	----

Warranty:

Frame	5 Years
Drive Train	13 Months
Electronics ²	13 Months
Batteries	12 Months

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.



High-Tech Styling Includes Highly Visible Side & Rear Lighting



Shown in Black Carbon Fiber Pattern



Power Folding



Right or Left Hand Controller



Side Lighting



Omni-Directional Wheels w/ Self-Cleaning V-Groove



Ergonomic Seating

Optional Accessory:



Wireless Remote



Phone: 800-265-9615
www.goldentechnologies.ca