



Venice PL447




Comfort Zone
by GOLDEN



Step into everyday luxury with the Venice by Golden. Designed for total relaxation, the Venice features a generously padded seam back, plush chaise pad, and a sleek open-arm design that includes a built-in cup holder, hidden storage, a power port with a removable reading light and tray table, a wireless phone charger, and a USB port—all within arm's reach. For added wellness, Nirvana™ Swedish Air Massage with Far Infrared heat comes standard.

Experience personalized comfort with our 5-Zone Comfort Zone™ system, enhanced with patented Twilight Tilt Technology and MaxiComfort Positioning. Effortlessly adjust the footrest, backrest, headrest, lumbar support, and tilt using our intuitive, programmable AutoDrive™ hand control.

Standard Fabrics:



Sand (ASD)



Stone (AST)

Standard with

Standard features:

- Patented Twilight™ Tilt Technology
- Patented MaxiComfort™ Positioning™ with Four-Way Motion Seat
- Adjustable Headrest & Lumbar
- Programmable AutoDrive™ Hand Control
- Nirvana™ Swedish Air Massage with Far Infrared Heat Therapy
- USB Charging Port
- Wireless Smartphone Charger
- Built-in Cup Holder
- Removable LED Reading Light & Tray Table
- Storage Compartments
- SmartTek® Diagnostic System
- Full & Straight Lift Options

PL447 Venice

Model Number	PL447
Size	Medium
Class II Medical Device	Yes
Weight Capacity	400 LB
Positioning	Infinite
Comfort Zone(s)	5
Twilight Tilt Technology	Yes
Chaise Seat	Yes
Pocketed Coil Seat	Yes
Back Type	Seam
Removable Back Type	Bracket
Overall Width	37"
Overall Height	45"
Overall Depth (Seated)	37"
Seat Width	22"
Floor to Top of Seat	20"
Seat Depth	21"
Seat to Top of Back	29"
Seat to Top of Arm	7"
Distance Required from Wall Reclined	25"
Back of Seat to End of Footrest Extended	43"
Lifetime Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism	Yes
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	Yes
Warranty on Mechanical Labor	1 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric, Foam and Pocket Coils	1 YR
Head Rest Cover	No
Arm Covers	No
Straight Lift	Yes
Hand Control Routing	Compartment
Extra Pocket Available	No
Nirvana Air Massage and Heat Available	Standard
Heat & Vibrating Massage Available	No
HeatWave Technology Available	No
Articulating Headrest Available	Standard
Adjustable Lumbar Available	Standard
Footrest Extension Available	No
USB Charging Port	Compartment
Battery Backup	Yes
Place of Origin	China

Standard Features & Accessories



Open Arm Design



LED Light



Cup Holder



Wireless Charger & Storage Compartment



USB Charging Port & Storage Compartment



Nirvana™ combines soothing Swedish massage with deep Far Infrared heat to relax and stretch muscles, relieving stiffness from prolonged sitting and restoring full-body wellness.

Smart Tek Diagnostic System STANDARD on all UltraComfort models

Comfort Zone™
by GOLDEN



5 Comfort Zones create the ultimate wellness comfort experience with our patented Twilight™ Tilt and MaxiComfort™ Positioning. Enjoy full-body support, including headrest and lumbar alignment, for effortless feet elevation and pressure-free comfort in every position.



WE PATENTED COMFORT™



Tel: 800-265-9615
www.goldentechnologies.ca

Designed and engineered in Old Forge, PA. Assembled in China.

Literature is current at time of printing. UltraComfort reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

Like us on Facebook
Facebook.com/GoldenTechCanada

LIT-SHT-PL447-CAN-0326