



Aurora PL480



Comfort Zone
by GOLDEN

The AURORA Power Lift Recliner combines timeless style with exceptional comfort. Featuring a supportive biscuit backrest and a plush, cushioned seat, it's designed for hours of relaxation. Available in three sizes, this recliner offers a seamless reclining experience with a simple two-button hand control, effortlessly transitioning you from sitting to reclining to standing—all while providing full-body support.

Standard Features:

- USB Charging Port
- SmartTek® Diagnostic System
- Available in Three Sizes: Small, Medium & Large
- Full, Mid & Straight Lift Options

Standard Fabrics



Anchor (EAC)



Sandstorm (ESS)

PL480 Aurora

Model Number	PL480	PL480	PL480
Size	Small	Medium	Large
Class II Medical Device	Yes	Yes	Yes
Weight Capacity	400 LB	400 LB	400 LB
Positioning	3-position	3-position	3-position
Comfort Zone(s)	1	1	1
Twilight Technology	No	No	No
Chaise Seat	Yes	Yes	Yes
Pocketed Coil Seat	No	No	No
Back Type	Biscuit	Biscuit	Biscuit
Removable Back Type	Bracket	Bracket	Bracket
Overall Width	31"	33"	35"
Overall Height	40.5"	42.5"	45.5"
Overall Depth (Seated)	37"	38"	41.5"
Seat Width	20"	21.5"	23.5"
Floor to Top of Seat	18"	19"	20.5"
Seat Depth	18.5"	20.5"	22"
Seat to Top of Back	27"	29"	31"
Seat to Top of Arm	7.5"	6.5"	7"
Distance Required from Wall Reclined	15"	17"	19"
Back of Seat to End of Footrest Reclined	41"	41"	42"
Lifetime Warranty, Chair Frame, Lift Frame, Recline Frame, Recline Mechanism	Yes	Yes	Yes
Warranty on Electrical Parts	3 YR	3 YR	3 YR
Pro Rated Years 4 Thru 7	Yes	Yes	Yes
Warranty on Mechanical Labor	1 YR	1 YR	1 YR
Warranty on Electrical Labor	3 YR	3 YR	3 YR
Warranty on Fabric, Foam, and Pocketed Coils	1 YR	1 YR	1 YR
Head Rest Cover	No	No	No
Arm Covers	No	No	No
Straight Lift	Yes	Yes	Yes
Hand Control Routing	Through the Pocket	Through the Pocket	Through the Pocket
Extra Pocket Available	No	No	No
Nirvana Air Massage & Heat Available	No	No	No
Heat Wave Technology Available	No	No	No
Heat & Vibrating Massage Available	No	No	No
Articulating Headrest Available	No	No	No
Adjustable Lumbar Available	No	No	No
Footrest Extension Available	No	No	No
USB Charging Port	Yes	Yes	Yes
Battery Backup	Yes	Yes	Yes
Place of Origin	China	China	China



Comfort Zone™
by GOLDEN

Comfort Zone 1 delivers a traditional recline experience with a simple two-button hand control that easily glides you from sitting to reclining to standing all while enjoying full body support.



WE PATENTED COMFORT™



Tel: 800-265-9615
goldentechnologies.ca

Like us on Facebook
Facebook.com/GoldenTechCanada

Designed and engineered in Old Forge, PA. Assembled in China.

Literature is current at time of printing. UltraComfort reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68

SHT-PL480-CAN-0926